

CLEANING GUIDE FOR HCL® Rx Fill Tray (#17487)

Materials:

Polystyrene

Recommended Disinfectant Options:

18627 Sterile Isopropanol Alcohol 32oz.

20142 Sterile HCL® Presaturated Polypropylene Wipes, 70% IPA, 9 x11

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7

19955 PREempt™ One-Step Surface Cleaner and Disinfectant, Trigger Spray, 32 oz.

19183 Sterile HYPO-CHLOR® 5.25% Trigger Spray, 16 oz.

20011 Sani-Cloth® Bleach Germicidal Wipes, 6 x 10½

19179 Sterile DECON-QUAT® SIMPLEMENTIX® Trigger Spray, 16 oz.

20010 Super Sani-Cloth® Germicidal Wipes, 6 x 7

18962 Sporidicin® Disinfectant Trigger Spray, 32 oz.

18963 Sporidicin® Disinfectant, 1-Gallon

Alternate Disinfectant Options:

3% Hydrogen Peroxide Based Cleaners

15% Sodium Hypochlorite Based Cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com

© Health Care Logistics®, Inc. 2025